

QA3 **DRIVE**™

◀800▶ Bariatric Stretcher



Specification Sheet

AneticAid® 

QA3 **DRIVE**™

◀800▶ Bariatric Stretcher

Extended platform width

Platform width has been extended to 800mm to suit bariatric patients.

Inbuilt motorised drive

No need to push or pull the stretcher: our inbuilt motorised DRIVE feature linked to the stretcher's fifth wheel offers superb manoeuvrability and control

Pushing handles

Four ergonomically positioned pushing handles for supreme manoeuvrability

Effortless travel

The QA3™ DRIVE <800> Bariatric Stretcher can travel up to gradients of 6°.

Long-life battery

Capable of **20 kilometers on a single charge**, the DRIVE feature is powered by an advanced nickel metal hydride rechargeable battery.

Simultaneous braking low rolling resistance castors

Polyurethane tyres offer low rolling resistance for high manoeuvrability: 150mm diameter central total lock castors fitted as standard, helping to offer a rigid and reliable work platform.

Exceptional patient comfort

Low height of 535mm (excluding mattress) & store away side rails allow for easier mount/dismount, particularly for frail and elderly patients. Also including a versatile positioning headrest, wide patient platform and pressure care mattress for comfort.

Press on / Press off 5th wheel steer

Our unique mechanism ensures maximum manoeuvrability and means the QA3 is easily guided by one support worker.

Fixed IV pole

Innovative quick release safety locking mechanism, integrated into the head end side rail for easy access and storage.

Pressure care mattress

The mattress is anti-microbial, water-resistant, vapour permeable, latex-free, and five per cent resistant to chlorine, ensuring optimal patient comfort and hygiene. An additional feature is a high-density foam layer underneath the main mattress structure, providing extra support.



SPECIFICATIONS		Code: 21115	
Minimum height	*535mm	Mattress length	2025mm
Maximum height	*835mm	Mattress width	800mm
*excludes mattress		Mattress depth	105mm
Minimum length	1700mm	Weight	167kg
Maximum length	2150mm	Backrest	0 to +90°
Overall width	905mm	Trendelenburg	0 to -12°
Castor diameter	150mm	Reverse Trendelenburg	*0 to -10°
Maximum patient weight		320kg	

*NOTE - Reverse Trendelenburg enabled through a build option

For more information on any of our products or service contracts, call:

sales@aneticaid.com **+44 (0)1943 878647** aneticaid.com



DESIGN & INNOVATION



MADE IN THE UK



GLOBAL DISTRIBUTION

Innovative **Medical** Technology - **Practically Applied**



Anetic Aid Ltd, Sapper Jordan Rossi Park, Otley Road, Baildon BD17 7AX, UK

AneticAid®